


3dk.berlin Shield Mask v4



Anti-fogging




Antibacterial



For your personal mask you need :
Scissors & Button Hole Elastic Band

I



II





Step 1

Insert the foil into the upper and lower part of the mask. The tabs on the top part will bend outwards and the bottom one's embed in the bottom part.


Step 2

Take the button hole elastic band and attach it to one side. Thread the band through the tabs and adjust the length to fit your head size. Then attach it to the other side.



Optional Step 3

You can increase the comfort, by sticking some foam rubber to the headband.



Shield **Mask**_{V4}

This is not a 100% shield against the CoVID19 virus, it only covers your mucous membranes from direct droplet infection.

It still applies:

keep your distance
wash your hands regularly
avoid gatherings

**Your 3dk.berlin team wishes everyone stamina and hopes for a speedy improvement.
www.3dk.berlin - masken@3dk.berlin - Ruppinerchaussee 325 13503 Berlin - 030 54909613**

ATTENTION: The Shield Mask V2 we supply have not been disinfected during production. If necessary, please disinfect them before first use.